



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08888 - Goat Anti-PCSK5 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PCSK5, proprotein convertase subtilisin/kexin type 5, PC5, PC6, PC6A, SPC6, OTTHUMP0000063532, prohormone convertase 5, proprotein convertase PC5, protease PC6, subtilisin/kexin-like protease PC5

Official Symbol: PCSK5

Accession Number(s): NP_006191.2

Human GeneID(s): [5125](#)

Immunogen

Peptide with sequence PDGSYQDTKKNLCKRK, from the internal region of the protein sequence according to NP_006191.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Preliminary experiments in lysates of cell lines K562 and MOLT4 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human