



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05467 - Goat Anti-PGLYRP Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PGLYRP, peptidoglycan recognition protein, PGRP, TAG7, PGRP-S, TNFSF3L, TNF superfamily, member 3 (LTB)-like (peptidoglycan recognition protein)

Official Symbol: PGLYRP

Accession Number(s): NP_005082

Human GeneID(s): [8993](#)

Immunogen

Peptide with sequence C-YHLIQNWPHYRSP, from the C Terminus of the protein sequence according to NP_005082.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 18-22kDa band seen in MOLT-4 lysate.

Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human