

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08421 - Goat Anti-Pleiotrophin Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PTN, pleiotrophin (heparin binding growth factor 8, neurite growth-promoting factor 1), HARP, HBGF8, HBNF, NEGF1, heparin affin regulatory protein, heparin binding growth factor 8, heparin-binding growth-associated molecule, neurite growth-promoting factor 1, pleiotrophin, HBBM, HB-GAM, HBGF-8, HBNF-1, heparin-binding brain mitogen, heparin-binding growth factor 8, heparin-binding neurite outgrowth-promoting factor 1, OSF-1, osteoblast-specific factor 1

Official Symbol: PTN

Accession Number(s): NP_002816.1

Human GeneID(s): 5764

Non-Human GenelD(s): 19242 (mouse), 24924 (rat)

Immunogen

Peptide with sequence C-KEGKKQEKMLD, from the C Terminus of the protein sequence according to NP 002816.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.


Western blot: Approx 18kDa band observed in lysates of Mouse fetal Brain (calculated MW of 18.9kDa according to Mouse NP_032999.1). Recommended concentration: 0.1-0.3µg/ml.

IHC: In paraffin embedded Human Prostate shows distinct staining at the base of the glands. Recommended concentration, 3-5µg/ml.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human, Mouse, Rat



EB08421 (3.8 μ g/ml) staining of paraffin embedded Human Prostate. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



EB08421 (0.1μg/ml) staining of Mouse fetal Brain (35μg protein in RIPA buffer). Primary incubation was 1 hour.

Detected by chemiluminescence.