



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05539 - Goat Anti-PMS1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PMS1, PMS1 postmeiotic segregation increased 1 (S. cerevisiae), PMSL1, hPMS1, HNPCC3, postmeiotic segregation 1, human homolog of yeast mutL, MLH1, DKFZp781M0253, FLJ98259, mismatch repair gene PMSL1, postmeiotic segregation 1, rhabdomyosarcoma antigen MU-RMS-40.10B, rhabdomyosarcoma antigen MU-RMS-40.10E

Official Symbol: PMS1

Accession Number(s): NP_000525.1; NP_001121615.1; NP_001121616.1

Human GeneID(s): [5378](#)

Important Comments: This antibody is expected to recognise all reported isoforms.

Immunogen

Peptide with sequence C-RPFFHHLTYLPETT, from the C Terminus of the protein sequence according to NP_000525.1; NP_001121615.1; NP_001121616.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Species Reactivity

Tested:

Expected from sequence similarity: Human