

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05539 - Goat Anti-PMS1 Antibody

Size: 100µg specific antibody in 200µl

Target Protein

Principal Names: PMS1, PMS1 postmeiotic segregation increased 1 (*S. cerevisiae*), PMSL1, hPMS1, HNPCC3, postmeiotic segregation 1, human homolog of yeast mutL, MLH1, DKFZp781M0253, FLJ98259, mismatch repair gene PMSL1, postmeiotic segregation 1, rhabdomyosarcoma antigen MU-RMS-40.10B, rhabdomyosarcoma antigen MU-RMS-40.10E

Official Symbol: PMS1

Accession Number(s): NP_000525.1; NP_001121615.1; NP_001121616.1

Human GeneID(s): [5378](#)

Important Comments: This antibody is expected to recognise all reported isoforms.

Immunogen

Peptide with sequence C-RPFFHHLTYLPETT, from the C Terminus of the protein sequence according to NP_000525.1; NP_001121615.1; NP_001121616.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: No signal obtained yet but low background observed in Human Brain and Kidney extracts at up to 2µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human