



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB11528 - Goat Anti-POFUT1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FUT12, GDP-fucose protein O-fucosyltransferase 1, KIAA0180, MGC2482, o-fucosyltransferase protein, O-Fuc-T, O-FucT-1, O-FUT, OTTHUMP00000030583, OTTHUMP00000030584, peptide-O-fucosyltransferase 1, POFUT1, protein O-fucosyltransferase 1

Official Symbol: POFUT1

Accession Number(s): NP_056167.1; NP_758436.1

Human GeneID(s): [23509](#)

Non-Human GeneID(s): 140484 (mouse), 311551 (rat)

Important Comments: This antibody is expected to recognize both reported isoforms (NP_056167.1; NP_758436.1).

Immunogen

Peptide with sequence C-TNLHVSQYQYFKLEP, from the internal region of the protein sequence according to NP_056167.1; NP_758436.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: In transfected HEK293 transiently expressing Human POFUT1 (myc and DYKDDDDK tagged), a band of approx. 48kDa is observed. No bands are observed in mock-transfected HEK293 and the same band is observed using anti-myc antibody. Recommended concentration, 0.3-1µg/ml.


IHC: In paraffin embedded Human Cerebral Cortex shows staining of the cytoplasm in neuronal bodies. Recommended concentration, 2-4µg/ml.

Purchase the +ve control! The transfected lysate used in QC is available for sale. Please email sales@everestbiotech.com to request datasheet, pricing and delivery information.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human POFUT1 with C-terminal MYC tag probed with EB11528 (1ug/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with EB11528 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.



EB11528 (2µg/ml) staining of paraffin embedded Human Cerebral Cortex. Steamed antigen retrieval with citrate buffer pH 6, HRP-staining.