



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10021 - Goat Anti-POU2F2 / OCT2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: homeobox protein, OCT2, Oct-2, OTF2, POU class 2 homeobox 2, POU domain, class 2, transcription factor 2, POU2F2

Official Symbol: POU2F2

Accession Number(s): NP_002689.1

Human GeneID(s): [5452](#)

Immunogen

Peptide with sequence KQGLDSPSEHTDTER-C, from the N Terminus of the protein sequence according to NP_002689.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 60kDa band observed in Human Tonsil and Mouse Thymus lysates (calculated MW of 49.5kDa according to NP_002689.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Pelis et al, Am J Physiol Renal Physiol. 2006 May;290(5):F1118-26. PMID: 16368738). Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig



EB10021 (0.1 μ g/ml) staining of Human Tonsil lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.