

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB09407 - Goat Anti-PPP2R4 / PP2A Antibody

Size: 100µg specific antibody in 200µl



#### **Target Protein**

**Principal Names:** PPP2R4, protein phosphatase 2A activator, regulatory subunit 4, MGC2184, PP2A, PR53, PTPA, OTTHUMP00000022333, PP2A phosphatase activator, PP2A, subunit B\1, phosphotyrosyl phosphatase activator, protein phosphatase 2A, regulatory subunit B\1, protein phosphatase 2A, regulatory subunit B\1 (PR 53)

Official Symbol: PTPA

Accession Number(s): NP\_821068.1; NP\_066954.2; NP\_821070.1

Human GeneID(s): 5524

Non-Human GenelD(s): 110854 (mouse), 362102 (rat)

**Important Comments:** This antibody is expected to recognize all reported isoforms (NP\_821068.1; NP\_066954.2; NP\_821070.1). Reported variants NP\_066954.2 and NP\_821067.1 represent identical protein.

## **Immunogen**

Peptide with sequence C-GSLLPIHPVTSG, from the C Terminus of the protein sequence according to NP\_821068.1; NP\_066954.2; NP\_821070.1.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

**Western blot:** Approx 37kDa band observed in Human, Mouse and Rat Liver lysates and in lysates of cell line HepG2 (calculated MW of 36.8kDa according to NP\_066954.2).


Recommended concentration: 0.5-2 $\mu$ g/ml. Primary incubation was 1 hour.

 $\textbf{IHC:} \ Paraffin \ embedded \ Human \ Testis. \ Recommended \ concentration: 3.75 \mu g/ml.$ 


## **Species Reactivity**

Tested: Human, Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB09407 (0.5 $\mu$ g/ml) staining of Human Liver lysate (lane 1), Mouse Liver (lane 2) and Rat Liver (lane 3) (35 $\mu$ g protein in RIPA buffer). Detected by chemiluminescence.



EB09407 (3.75 $\mu$ g/ml) staining of paraffin embedded Human Testis. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.