



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12266 - Goat Anti-prealbumin / transthyretin (aa62-74) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TTR, transthyretin, CTS, CTS1, HsT2651, PALB, TBPA, ATTR, carpal tunnel syndrome 1, prealbumin, amyloidosis type I, thyroxine-binding prealbumin

Official Symbol: TTR

Accession Number(s): NP_000362.1

Human GeneID(s): [7276](#)

Immunogen

Peptide with sequence C-EPFASGKTSESGE, from the internal region of the protein sequence according to NP_000362.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 16kDa band observed in Human Albumin-depleted Serum (calculated MW of 15.9kDa according to NP_000362.1). Recommended concentration: 0.01-0.03µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Cow



EB12266 (0.03 μ g/ml) staining of Human Albumin-depleted Serum (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.