

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07025 - Goat Anti-PTCH (C Terminus) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PTCH, patched homolog (Drosophila), HGNC:9585, BCNS, HPE7, NBCCS, PTC, PTC1, patched, patched (Drosophila) homolog, PTCH1, patched homolog 1 (Drosophila), RP11-435O5.3, FLJ26746, FLJ42602, PTCH11, OTTHUMP00000021709, OTTHUMP00000021710, PTCH protein +12b, PTCH protein +4', PTCH protein -10

Official Symbol: PTCH1

Accession Number(s): NP_000255.1; NP_001077072.1; NP_001077071.1; NP_001077073.1; NP_001077074.1; NP_001077075.1; NP_001077076.1

Human GeneID(s): [5727](#)

Important Comments: This antibody is expected to recognise all reported isoforms. Variants (NP_001077073.1; NP_001077074.1; NP_001077075.1; NP_001077076.1) encode the same isoform.

Immunogen

Peptide with sequence ELQDVECEERPR, from the C Terminus of the protein sequence according to NP_000255.1; NP_001077072.1; NP_001077071.1; NP_001077073.1; NP_001077074.1; NP_001077075.1; NP_001077076.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

IHC: Paraffin embedded Human Brain (Coretx). Recommended concentration: 8µg/ml.

Immunofluorescence: Strong expression of the protein seen in the Golgi of HeLa and NIH3T3 cells. Recommended concentration: 10µg/ml.

Species Reactivity

Tested: Human, Mouse

Expected from sequence similarity: Human, Mouse, Rat

Specific Reference


This antibody (previous batch) has been successfully used in IHC and IF on Human:

Abiria SA, Williams TV, Munden AL, Grover VK, Wallace A, Lundberg CJ, Valadez JG, Cooper MK.


Expression of Hedgehog ligand and signal transduction components in mutually distinct isocitrate dehydrogenase mutant glioma cells supports a role for paracrine signaling.

J Neurooncol. 2014 May 28.


PMID: 24867209




EB07025 Immunofluorescence analysis of paraformaldehyde fixed HeLa cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing Golgi/cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).



EB07025 Immunofluorescence analysis of paraformaldehyde fixed NIH3T3 cells, permeabilized with 0.15% Triton. Primary incubation 1 hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing Golgi staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).



EB07025 (8µg/ml) staining of paraffin embedded Human Cortex. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



EB07025 Negative Control showing staining of paraffin embedded Human Cortex, with no primary antibody.