



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05262 - Goat Anti-PUMA Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BBC3, BCL2 binding component 3, JFY1, PUMA, PUMA/JFY1, Bcl-2 binding component 3, Bcl-2 binding component 3 PUMA/JFY1 protein

Official Symbol: BBC3

Accession Number(s): NP_055232.1

Human GeneID(s): [27113](#)

Immunogen

Peptide with sequence C-LPRGHRAPEMEPN, from the C Terminus of the protein sequence according to NP_055232.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Western blot: Approx 22kDa band seen in NIH-3T3 lysates (Predicted MW of 21kDa according to NP_055232). Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human