

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05262 - Goat Anti-PUMA Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BBC3, BCL2 binding component 3, JFY1, PUMA, PUMA/JFY1, Bcl-2

binding component 3, Bcl-2 binding component 3 PUMA/JFY1 protein

Official Symbol: BBC3

Accession Number(s): NP_055232.1

Human GeneID(s): 27113

Immunogen

Peptide with sequence C-LPRGHRAPEMEPN, from the C Terminus of the protein sequence according to NP_055232.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Western blot: Approx 22kDa band seen in NIH-3T3 lysates (Predicted MW of 21kDa according to NP_055232). Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human