



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05262 - Goat Anti-PUMA Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BBC3, BCL2 binding component 3, JFY1, PUMA, PUMA/JFY1, Bcl-2 binding component 3, Bcl-2 binding component 3 PUMA/JFY1 protein

Official Symbol: BBC3

Accession Number(s): NP_055232.1

Human GeneID(s): [27113](#)

Immunogen

Peptide with sequence C-LPRGHRAPEMEPN, from the C Terminus of the protein sequence according to NP_055232.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Western blot: Approx 22kDa band seen in NIH-3T3 lysates (Predicted MW of 21kDa according to NP_055232). Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human