



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB2ND-51-3001 - Goat anti-Rabbit IgG antibody

Size: 500µg



Immunogen

Purified rabbit whole IgG (unconjugated).

Purification

Ammonium sulphate precipitation of serum followed by antigen affinity purification.

Buffer

Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative 10 mM Tris, PH 8.2, with 150 mM NaCl, 10 mg/ml BSA and 0.1% sodium azide. Quantity: 0.2mg in 0.5mg/ml.

Storage and Handling

Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles.

Applications

Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. Recommended dilutions:

ELISA: 1:128000 WB: 1:3000 to 1:10000. Please optimize dilutions for your particular assay conditions. Starting test dilution: 0.5µg/ml for ELISA and WB.