



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06053 - Goat Anti-Rabenosyn 5 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ZFYVE20, Rabenosyn 5, FYVE-finger-containing Rab5 effector protein rabenosyn-5, zinc finger, FYVE domain containing 20, zinc finger, FYVE domain containing 20, FLJ34993, MGC126210, Rabenosyn-5

Official Symbol: RBSN

Accession Number(s): NP_071735.2

Human GeneID(s): [64145](#)

Immunogen

Peptide with sequence C-ELKHTLAKQKGGTD, from the C Terminus of the protein sequence according to NP_071735.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Approx 110 kDa band observed in A431, Fibroblast (epidermoid carcinoma) lysates (calculated MW of 88.8 kDa according to NP_071735). Band migration is consistent with literature (PMID 11062261, 15020713). Recommended concentration: 0.3 - 1.5 µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB06053 (0.3 μ g/ml) staining of A431 lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.