

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07140 - Goat Anti-Rabring 7 / BCA2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Rabring 7, BCA2, ZNF364, zinc finger protein 364, CL469780, RNF115, breast cancer-associated gene 2, ring finger protein 115

Official Symbol: RNF115

Accession Number(s): NP_055270.1

Human GeneID(s): [27246](#)

Non-Human GeneID(s): 67845 (mouse), 362002 (rat)

Immunogen

Peptide with sequence C-DQDNRRANERGHQTH, from the internal region of the protein sequence according to NP_055270.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 35kDa band observed in Human lysates of cell line HEPG2 (calculated MW of 33.7kDa according to NP_055270.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog



EB07140 (0.3 μ g/ml) staining of HEPG2 cell lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.