

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06371 - Goat Anti-RASGRP3 / GRP3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GRP3, KIAA0846, guanine nucleotide exchange factor for Rap1, RAS

guanyl releasing protein 3 (calcium and DAG-regulated)

Official Symbol: RASGRP3

Accession Number(s): NP_733772.1

Human GeneID(s): 25780

Immunogen

Peptide with sequence GSSGLGKAATLDE-C, from the N Terminus of the protein sequence according to NP_733772.1.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Western Blot: Preliminary experiments gave an approx 48kDa band in Human Brain and Kidney lysates at 1ug/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 79kDda according to NP_733772.1. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human