

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB10077 - Goat Anti-RDH5 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: 11-cis RDH, 9-cis-retinol specific dehydrogenase, FLJ39337, FLJ97089, HSD17B9, RDH1, retinol dehydrogenase 1, retinol dehydrogenase 5 (11-cis

and 9-cis), retinol dehydrogenase 5 (11-cis/9-cis), SDR9C5, short chain

dehydrogenase/reductase family 9C,, RDH5

Official Symbol: RDH5

Accession Number(s): NP_002896.2

Human GeneID(s): 5959

Non-Human GenelD(s): 366791 (rat)

Immunogen

Peptide with sequence C-RTPVTNLESLEKT, from the internal region of the protein sequence according to NP_002896.2.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Rat, Cow