

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10077 - Goat Anti-RDH5 Antibody

Size: 100µg specific antibody in 200µl

Target Protein

Principal Names: 11-cis RDH, 9-cis-retinol specific dehydrogenase, FLJ39337, FLJ97089, HSD17B9, RDH1, retinol dehydrogenase 1, retinol dehydrogenase 5 (11-cis and 9-cis), retinol dehydrogenase 5 (11-cis/9-cis), SDR9C5, short chain dehydrogenase/reductase family 9C., RDH5

Official Symbol: RDH5

Accession Number(s): NP_002896.2

Human GeneID(s): [5959](#)

Non-Human GeneID(s): 366791 (rat)

Immunogen

Peptide with sequence C-RTPVTNLESLEKT, from the internal region of the protein sequence according to NP_002896.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave bands at approx 70kDa and 40kDa in Mouse Eye lysates after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the bands we observe given the calculated size of 35.0kDa according to NP_002896.2. Both detected bands were successfully blocked by incubation with the immunizing peptide (and BLAST results with the immunizing peptide sequence did not identify any other proteins to explain the additional bands).

Species Reactivity

Tested:

Expected from sequence similarity: Human, Rat, Cow