

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB05219 - Goat Anti-RIN1 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

Principal Names: RIN1, Ras and Rab interactor 1, ras inhibitor 1

Official Symbol: RIN1

Accession Number(s): NP\_004283

Human GenelD(s): 9610

### **Immunogen**

Peptide with sequence C-PGEPEAEGSRAAEE, from the C Terminus of the protein sequence according to NP\_004283.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

# **Applications Tested**

**Western blot:** Western blot: Approx 100kDa band and 35kDa bands seen in HeLa cell lysates. Recommended for use at 1-3µg/ml.

Please note that the stronger 100 kDa band corresponds to the reference sequence NP\_004283, whereas the fainter 35 kDa band may correspond to a novel isoform of this protein.

## **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human