

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05390 - Goat Anti-RIPK2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: RIPK2, receptor-interacting serine-threonine kinase 2, RICK, RIP2,

CARD3, CARDIAK, receptor interacting protein 2

Official Symbol: RIPK2

Accession Number(s): NP_003812

Human GeneID(s): 8767

Immunogen

Peptide with sequence C-SRSPSLNLLQNKSM, from the C Terminus of the protein sequence according to NP_003812.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 60-65kDa band seen in Daudi lysate. Recommended for use at 1-3μg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human