



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05004 - Goat Anti-RIPK3 / RIP3 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** RIPK3, receptor-interacting serine-threonine kinase 3, RIP3, receptor-interacting protein 3, RIP-like protein kinase 3

**Official Symbol:** RIPK3

**Accession Number(s):** NP\_006862.2

**Human GeneID(s):** [11035](#)

### Immunogen

Peptide with sequence C-AWSRPQGWNHSGK, from the C Terminus of the protein sequence according to NP\_006862.2.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Western blot:** Western blot: Approx 45-50kDa band seen in HeLa, Jurkat and NSO (mouse) cell lysates. Recommended for use at 1-3 µg/ml. Overnight incubation is recommended with this antibody.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human