



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07774 - Goat Anti-Robo3 (mouse) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Robo3, roundabout homolog 3 (Drosophila), Rbig1, RIG1, Rig-1, Robo3a, Robo3b, retinoblastoma inhibiting gene 1

Official Symbol: Robo3

Accession Number(s): NP_035378.1; AAD11628.1

Non-Human GeneID(s): 19649 (mouse), 315564 (rat)

Important Comments: Both variants represent identical protein

Immunogen

Peptide with sequence C-RRQKQRKELSHYT, from the internal region of the protein sequence according to NP_035378.1; AAD11628.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave an approx 75-80kDa band in Mouse Brain and Rat Brain lysates after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 143kDa according to NP_035378.1 and AAD11628.1. The 75-80kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Mouse, Rat