



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09202 - Goat Anti-SCGB3A1 / UGRP2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SCGB3A1, secretoglobin, family 3A, member 1, HIN-1, HIN1, LU105, MGC87867, PnSP-2, UGRP2, cytokine high in normal-1, pneumo secretory protein 2

Official Symbol: SCGB3A1

Accession Number(s): NP_443095.2

Human GeneID(s): [92304](#)

Immunogen

Peptide with sequence NHLIEGSQKCVAE, from the internal of the protein sequence according to NP_443095.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave an approx 33kDa band in Human Breast cancer and normal Lung lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 10.1kDa according to NP_443095.2. The 33kDa band was successfully blocked by incubation with the immunizing peptide.

Species Reactivity

Tested:

Expected from sequence similarity: Human