

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB09202 - Goat Anti-SCGB3A1 / UGRP2 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** SCGB3A1, secretoglobin, family 3A, member 1, HIN-1, HIN1, LU105, MGC87867, PnSP-2, UGRP2, cytokine high in normal-1, pneumo secretory protein 2

Official Symbol: SCGB3A1

Accession Number(s): NP\_443095.2

Human GeneID(s): 92304

### **Immunogen**

Peptide with sequence NHLIEGSQKCVAE, from the internal of the protein sequence according to NP\_443095.2.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

**Western blot:** Preliminary experiments gave an approx 33kDa band in Human Breast cancer and normal Lung lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 10.1kDa according to NP\_443095.2. The 33kDa band was successfully blocked by incubation with the immunizing peptide.

## **Species Reactivity**

Tested:

Expected from sequence similarity: Human