

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB10428 - Goat Anti-SEPT4 (aa33-44) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ARTS, BRADEION, CE5B3, cell division control-related protein 2, cerebral protein 7, H5, hCDCREL-2, hucep-7, MART, peanut-like 2, PNUTL2, SEP4,

septin 4

Official Symbol: SEPTIN4

Accession Number(s): NP_004565.1; NP_536340.1; NP_536341.1; NP_001185642.1

Human GeneID(s): 5414

Non-Human GenelD(s): 18952 (mouse)

Important Comments: This antibody is expected to recognize all reported isoforms (NP_004565.1; NP_536340.1; NP_536341.1, NP_001185642.1). Amino acid numbering

in name refers to NP_004565.1 sequence.

Immunogen

Peptide with sequence C-ELSKFVKDFSGN, from the internal region of the protein sequence according to NP_004565.1; NP_536340.1; NP_536341.1; NP_001185642.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.


Western blot: Approx 55kDa band observed in nuclear lysates of cell lines Jurkat and in LNCaP cell lysates, 50-55kDa in Human Cerebellum lysates, and approx. 50kDa in Mouse Brain lysates (calculated MW of 53.0kDa according to Human NP_536341.1 and 49.4kDa according to Mouse NP_001271323.1). Recommended concentration: 0.1-0.3μg/ml. Primary incubation 1 hour at room temperature.

Flow Cytometry: Flow cytometric analysis of Jurkat cells. Recommended concentration: 10ug/ml.

Species Reactivity


Tested: Human, Mouse

Expected from sequence similarity: Human, Mouse




EB10428 (0.3μg/ml) staining of Jurkat nuclear (A) and LNCaP (B) cell lysate (35μg protein in RIPA buffer).

Detected by chemilluminescence.



EB10428 (0.3μg/ml) staining of 1Human Cerebellum (A) and Mouse Brain (B) lysate (35μg protein in RIPA buffer). Detected by chemiluminescence.



EB10428 Flow cytometric analysis of paraformaldehyde fixed Jurkat cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control:

Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.