

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08539 - Goat Anti-SERPINB11 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SERPINB11, serpin peptidase inhibitor, clade B (ovalbumin), member 11, EPIPIN, MGC163242, MGC163244, serine (or cysteine) proteinase inhibitor, clade B

(ovalbumin), member 11, serpin peptidase inhibitor, clade B, member 11

Official Symbol: SERPINB11

Accession Number(s): NP_536723.2

Human GeneID(s): 89778

Immunogen

Peptide with sequence EFSQINQPDSNCT, from the internal region of the protein sequence according to NP_536723.2.

Please note the <u>peptide</u> is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Preliminary experiments in Human Lung, Placenta, Prostate and Tonsil lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human