

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08577 - Goat Anti-SERPINB6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SERPINB6, serpin peptidase inhibitor, clade B (ovalbumin), member 6, CAP, DKFZp686I04222, MGC111370, MSTP057, PI6, PTI, protease inhibitor 6 (placental thrombin inhibitor), serine (or cysteine) proteinase inhibitor, clade B (ovalbumin), member

6

Official Symbol: SERPINB6

Accession Number(s): NP_004559.4

Human GenelD(s): 5269

Immunogen

Peptide with sequence C-NWDEQFDKENTEER, from the internal region of the protein sequence according to NP_004559.4.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 40kDa band observed in Human Ovary lysates (calculated MW of 42.6kDa according to NP_004559.4). Recommended concentration: 0.1-.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB08577 (0.1 μ g/ml) staining of Human Ovary lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.