

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05292 - Goat Anti-SIAH1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SIAH1, seven in absentia homolog 1 (Drosophila), Siah-1, hSIAH1, HUMSIAH, Siah-1a, FLJ08065, seven in absentia homolog 1, siah E3 ubiquitin protein

ligase 1, SIAH1A

Official Symbol: SIAH1

Accession Number(s): NP_003022.3; NP_001006611.1

Human GeneID(s): 6477

Non-Human GenelD(s): 20437 (mouse), 140941 (rat)

Important Comments: This product is expected to recognize both reported isoforms

(NP_003022.3; NP_001006611.1)

Immunogen

Peptide with sequence SRQTATALPTGTSKC, from the N Terminus of the protein sequence according to NP_003022.3; NP_001006611.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx. 37kDa band observed in Human Liver lysates and approx. 30+37kDa bands in Rat Liver lysates (calculated MW of 34.6kDa according to Human NP_001006611.1 and 31.1kDa according to Rat NP_543181.2 and 38.2kDa according to Rat XP_038953364.1). Recommended concentration: 1-3ug/ml. Primary incubation 1 hour at room temperature.

Immunofluorescence: Strong expression of the protein seen in the nucleus and cytoplasm of U2OS and HepG2 cells. Recommended concentration: 10µg/ml.

Species Reactivity

Tested: Human, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow

Specific References

This antibody (previous batch) has been successfully used in WB and IP on Human:

Liman J, Sroka K, Dohm CP, Deeg S, Bähr M, Kermer P.

Modulation of Huntingtin Toxicity by BAG1 is Dependent on an Intact BAG Domain.

Molecules 2010, 15, 6678-6687

PMID: 0


This antibody (previous batch) has been successfully used in Western blot on Drosophila:

Sroka K, Voigt A, Deeg S, Reed JC, Schulz JB, Bähr M, Kermer P.

BAG1 modulates huntingtin toxicity, aggregation, degradation, and subcellular distribution.


J Neurochem. 2009 Nov;111(3):801-7.

PMID: 19712056




EB05292 optimised QC. Primary incubation 1 hour at room temperature.

Images A+B: Human Liver lysate + peptide, incubation at primary Ab concentration 1ug/ml, C+D: Rat Liver lysate + peptide, incubation at primary Ab concentration 3ug/ml. (Loaded 35µg protein in RIPA buffer, per lane). Detected by chemiluminescence.



EB05292 Immunofluorescence analysis of paraformaldehyde fixed U2OS cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing nuclear and cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody.



EB05292 Immunofluorescence analysis of paraformaldehyde fixed HepG2 cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing nuclear and cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody