

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB05279 - Goat Anti-SLAP2 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

Principal Names: SLA2, Src-like-adaptor 2, SLAP2, SLAP-2, FLJ21992, C20orf156,

dJ977B1.1, src-like-adaptor 2, Src-like adapter protein-2

Official Symbol: SLA2

Accession Number(s): NP\_115590; NP\_778252

Human GeneID(s): 84174

## **Immunogen**

Peptide with sequence C-ISLNDEAVSLDDA, from the C Terminus of the protein sequence according to NP\_115590; NP\_778252.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### **Applications Tested**

**Western blot:** Western blot: Approx 35kDa band seen in MOLT-4 cell lysate. Recommended for use at 2-4µg/ml.

## **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human