



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05130 - Goat Anti-SLC16A1 / MCT1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SLC16A1, solute carrier family 16 (monocarboxylic acid transporters), member 1, MCT, MCT1, Solute carrier family 16 (monocarboxylic acid transporters), solute carrier family 16, member 1 (monocarboxylic acid transporter 1), FLJ36745, HHF7, MGC44475, monocarboxylate transporter 1, solute carrier family 16, member 1

Official Symbol: SLC16A1

Accession Number(s): NP_003042.3

Human GeneID(s): [6566](#)

Immunogen

Peptide with sequence C-DTEGGPKKEESPV, from the C Terminus of the protein sequence according to NP_003042.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Western Blot: Approx 60kDa band observed in lysates of Human Burkitt's Lymphoma cell line DAUDI (calculated MW of 54.0kDa according to NP_003042.2).
Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human