

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

## EB05130 - Goat Anti-SLC16A1 / MCT1 Antibody

Size: 100µg specific antibody in 200µl



#### **Target Protein**

**Principal Names:** SLC16A1, solute carrier family 16 (monocarboxylic acid transporters), member 1, MCT, MCT1, Solute carrier family 16 (monocarboxylic acid transporters), solute carrier family 16, member 1 (monocarboxylic acid transporter 1), FLJ36745, HHF7, MGC44475, monocarboxylate transporter 1, solute carrier family 16, member 1

Official Symbol: SLC16A1

Accession Number(s): NP\_003042.3

Human GeneID(s): 6566

## **Immunogen**

Peptide with sequence C-DTEGGPKEEESPV, from the C Terminus of the protein sequence according to NP\_003042.3.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:1000.

**Western blot:** Western Blot: Approx 60kDa band observed in lysates of Human Burkitt's Lymphoma cell line DAUDI (calculated MW of 54.0kDa according to NP\_003042.2). Recommended concentration: 0.1-0.3µg/ml.

# **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human