



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05133 - Goat Anti-SLC16A4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SLC16A4, solute carrier family 16 (monocarboxylic acid transporters), member 4, MCT4, monocarboxylate transporter 4

Official Symbol: SLC16A4

Accession Number(s): NP_004687

Human GeneID(s): [9122](#)

Immunogen

Peptide with sequence C-FVPLAERWKNSLT, from the C Terminus of the protein sequence according to NP_004687.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 55-60kDa band seen in 293 cell lysate.

Recommended for use at 0.5-2 µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human