

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05133 - Goat Anti-SLC16A4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SLC16A4, solute carrier family 16 (monocarboxylic acid transporters),

member 4, MCT4, monocarboxylate transporter 4

Official Symbol: SLC16A4

Accession Number(s): NP_004687

Human GeneID(s): 9122

Immunogen

Peptide with sequence C-FVPLAERWKNSLT, from the C Terminus of the protein sequence according to NP_004687.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 55-60kDa band seen in 293 cell lysate. Recommended for use at $0.5-2 \mu g/ml$.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human