

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

#### www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB09132 - Goat Anti-SLC26A5 / prestin Antibody

Size: 100µg specific antibody in 200µl



### **Target Protein**

**Principal Names:** SLC26A5, solute carrier family 26, member 5 (prestin), DFNB61, MGC118886, MGC118887, MGC118888, MGC118889, PRES, OTTHUMP00000195086, deafness, neurosensory, autosomal recessive, 61, prestin, prestin (motor protein)

Official Symbol: SLC26A5

Accession Number(s): NP\_945350.1; NP\_996766.1; NP\_996767.1; NP\_996768.1

Human GeneID(s): 375611

Important Comments: This antibody is expected to recognize reported isoforms

NP\_996768.1; NP\_996767.1; NP\_996766.1; NP\_945350.1.

### **Immunogen**

Peptide with sequence ERLHTKDKVPDSIAD -C, from the N Terminus (near) of the protein sequence according to NP\_945350.1; NP\_996766.1; NP\_996767.1; NP\_996768.1.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

**Western blot:** Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.

#### **Species Reactivity**

Tested:

Expected from sequence similarity: Human