

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09132 - Goat Anti-SLC26A5 / prestin Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SLC26A5, solute carrier family 26, member 5 (prestin), DFNB61, MGC118886, MGC118887, MGC118888, MGC118889, PRES, OTTHUMP00000195086, deafness, neurosensory, autosomal recessive, 61, prestin, prestin (motor protein)

Official Symbol: SLC26A5

Accession Number(s): NP_945350.1; NP_996766.1; NP_996767.1; NP_996768.1

Human GeneID(s): [375611](#)

Important Comments: This antibody is expected to recognize reported isoforms NP_996768.1; NP_996767.1; NP_996766.1; NP_945350.1.

Immunogen

Peptide with sequence ERLHTKDKVPDSIAD -C, from the N Terminus (near) of the protein sequence according to NP_945350.1; NP_996766.1; NP_996767.1; NP_996768.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.

Species Reactivity

Tested:

Expected from sequence similarity: Human