

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05081 - Goat Anti-SLI Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SLI, neuronal Shc adaptor homolog, SCK, SHCB

Official Symbol: SHC2 Human GeneID(s): 25759

Immunogen

Peptide with sequence C-ESELHLRGVVSREP, from the C Terminus of the protein sequence according to .

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 60 kDa and 30ka bands seen in Jurkat Lysate. Recommended for use at 2-4 µg/ml.

Please note that we also see an approximately 30kDa band in addition to the expected band. This product may recognise a Shc-related protein, whose predicted MW is 28kDa, according to AAC97988.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human