



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07047 - Goat Anti-SLITRK1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SLITRK1, SLIT and NTRK-like family, member 1, HGNC:20297, KIAA0918, KIAA1910, LRRC12, RP11-395N17.1, leucine rich repeat containing 12, slit and trk like 1 protein, slit and trk like gene 1, FLJ54428

Official Symbol: SLITRK1

Accession Number(s): NP_443142.1

Human GeneID(s): [114798](#)

Immunogen

Peptide with sequence NADGAHRVYDCGSHS, from the C Terminus of the protein sequence according to NP_443142.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary experiments in human brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow