



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05149 - Goat Anti-SMAD4 / MADH4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MADH4, SMAD4, JIP, DPC4, MAD, mothers against decapentaplegic homolog 4 (Drosophila), Mothers against decapentaplegic, Drosophila, homolog of, 4, MAD (mothers against decapentaplegic, Drosophila) homolog 4, SMAD family member 4, OTTHUMP00000163548, SMAD, mothers against DPP homolog 4, deleted in pancreatic carcinoma locus 4, mothers against decapentaplegic homolog 4

Official Symbol: SMAD4

Accession Number(s): NP_005350.1

Human GeneID(s): [4089](#)

Non-Human GeneID(s): 50554 (rat)

Immunogen

Peptide with sequence C-HTMPIADPQPLD, from the C Terminus of the protein sequence according to NP_005350.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave a band at approx 75kDa in lysates of cell line NIH3T3. Please note that currently we cannot find an explanation in the literature

Immunoprecipitation: This antibody was deemed fit for IP under native conditions (observations from anonymous customer).

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig Cow