

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05150 - Goat Anti-SMAD6 / MADH6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MADH6, MAD, mothers against decapentaplegic homolog 6 (Drosophila), MADH7, SMAD6, Mothers against decapentaplegic, drosophila, homolog of,

6, MAD (mothers against decapentaplegic, Drosophila) homolog 6

Official Symbol: SMAD6 Human GenelD(s): 4091

Immunogen

Peptide with sequence PCWLEILLNNPR, from the C Terminus of the protein sequence according to .

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 55kDa band seen in Jurkat lysate. Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human