

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05268 - Goat Anti-SMURF1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SMURF1, E3 ubiquitin ligase SMURF1

Official Symbol: SMURF1

Accession Number(s): AAF08298

Human GeneID(s): 57154

Immunogen

Peptide with sequence DVRGSLQTPQNRPHC, from the internal region of the protein sequence according to AAF08298.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: No signal obtained yet but low background observed in HeLa and NIH-3T3 cell lysates at upto 2μg/ml.

Species Reactivity

Tested:

Expected from sequence similarity: Human