

#### **International Office**

**Everest Biotech Ltd** 

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB05675 - Goat Anti-SNAPAP Antibody

Size: 100µg specific antibody in 200µl



**Target Protein** 

Principal Names: SNAPAP, SNARE associated protein snapin, SNAPIN

Official Symbol: SNAPAP

Accession Number(s): NP\_036569

Human GeneID(s): 23557

### **Immunogen**

Peptide with sequence C-LDSGIYPPGSPGK, from the C Terminus of the protein sequence according to NP\_036569.

Please note the peptide is available for sale.

# Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

**Western blot:** Western blot: Approx 16kDa band seen in Human Brain lysate [Predicted MW of 15kDa according to NP\_036569]. Recommended for use at 1-3µg/ml.

#### **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human