

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB07042 - Goat Anti-SNM1 / DCLRE1A (internal) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: DCLRE1A, DNA cross-link repair 1A (PSO2 homolog, S. cerevisiae), HGNC:17660, KIAA0086, PSO2, SNM1, DNA-crosslink repair gene SNM1, SNM1A, DNA

cross-link repair 1A, OTTHUMP00000020516

Official Symbol: DCLRE1A

Accession Number(s): NP_055696.3

Human GeneID(s): 9937

Immunogen

Peptide with sequence PSPNQKRSSQCKRK, from the internal region of the protein sequence according to NP_055696.3.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.


Western blot: P

IHC: Paraffin embedded Human Skeletal Muscle. Recommended concentration: 5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



 $EB07042~(5\mu g/ml)~staining~of~paraffin~embedded~Human~Skeletal~Muscle.~Steamed~antigen~retrieval~with~citrate~buffer~pH~6,~AP-staining.$