

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06095 - Goat Anti-SNT2 / FRS3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: fibroblast growth factor receptor substrate 3, SNT-2, suc1-associated neurotrophic factor target 2 (FGFR signalling adaptor)

Official Symbol: FRS3

Accession Number(s): NP_006644

Human GeneID(s): 10817

Immunogen

Peptide with sequence GSCCSCLNRDSVPD, from the N Terminus of the protein sequence according to NP_006644.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Western Blot: Preliminary experiments gave no signal but low background in Jurkat or Hela lysates at up to $1\mu g/ml$. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human