

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB11617 - Goat Anti-SOX11 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

Principal Names: Sox11, SRY (sex determining region Y)-box 11

Official Symbol: Sox11

Accession Number(s): NP\_003099.1

Human GeneID(s): 6664

Non-Human GeneID(s): 20666 (mouse), 84046 (rat)

## **Immunogen**

Peptide with sequence C-DWLEANFSDLVFTY, from the C Terminus of the protein sequence according to NP\_003099.1.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:8000.

## **Species Reactivity**

Tested:

**Expected from sequence similarity:** Human, Mouse, Rat, Pig, Cow, Xenopus, Danio

rerio