

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09486 - Goat Anti-SP110 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SP110, SP110 nuclear body protein, FLJ22835, IFI41, IFI75, IPR1, VODI, interferon-induced protein 41, 30kD, interferon-induced protein 75, 52kD, phosphoprotein 41, phosphoprotein 75, transcriptional coactivator Sp110

Official Symbol: SP110

Accession Number(s): NP_004500.3; NP_004501.3; NP_536349.2

Human GeneID(s): [3431](#)

Important Comments: This antibody is expected to recognize all three reported isoforms (NP_004500.3; NP_004501.3; NP_536349.2).

Immunogen

Peptide with sequence C-PQEVSSTPSDKKGGK, from the internal region of the protein sequence according to NP_004500.3; NP_004501.3; NP_536349.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.


IHC: In paraffin embedded Human Tonsil shows nucleolar staining of lymphocytes.

Recommended concentration, 3-6µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB09486 (3.8µg/ml) staining of paraffin embedded Human Tonsil. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.