



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07160 - Goat Anti-STELLAR Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** STELLAR, germ and embryonic stem cell enriched protein STELLA, developmental pluripotency associated 3 pseudogene, stella-like, STELLA

**Official Symbol:** DPPA3

**Accession Number(s):** NP\_954980.1

**Human GeneID(s):** [400206](#)

### Immunogen

Peptide with sequence C-NKEPKGVKKESRP, from the internal region of the protein sequence according to NP\_954980.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:64000.

**Western blot:** Preliminary experiments gave an approx 26kDa band in human ovary and placenta lysates after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 17.9kDa according to NP\_954980.1. The 26kDa band was successfully blocked by incubation with the immunizing peptide.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human