



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB09334 - Goat Anti-STK36 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** STK36, serine/threonine kinase 36, fused homolog (Drosophila), DKFZp434N0223, FU, KIAA1278, fused, serine/threonine kinase 36

**Official Symbol:** STK36

**Accession Number(s):** NP\_056505.1

**Human GeneID(s):** [27148](#)

### Immunogen

Peptide with sequence C-QSLPHSSPRPASAK, from the C Terminus of the protein sequence according to NP\_056505.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:16000.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Dog