



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05596 - Goat Anti-TAF4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TAF4: TAF4 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 135 kD, TAF2C, TAF4A, TAF2C1, TAFII130, TAFII135, transcription initiation factor TFIID 135 kD subunit, TATA box binding protein (TBP)-associated factor, RNA polymerase II, C1, 130kD, TAF4A RNA polymerase II, TATA box binding protein (TBP)-associated factor, 135 kD

Official Symbol: TAF4

Accession Number(s): NP_003176

Human GeneID(s): [6874](#)

Non-Human GeneID(s): 117169 (mouse)

Immunogen

Peptide with sequence C-TSHSLLLYKAFLK, from the C Terminus of the protein sequence according to NP_003176.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 125kDa band seen in Human Brain lysate (Predicted MW of 119kDa according to NP_003176). Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse