

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB11633 - Goat Anti-TAS2R4 (aa157-171) Antibody

Size: 100µg specific antibody in 200µl

## **Target Protein**

Principal Names: candidate taste receptor T2R4, MGC163311,

OTTHUMP00000212311, T2R4, TAS2R4, taste receptor type 2 member 4, taste receptor,

type 2, member 4

Official Symbol: TAS2R4

Accession Number(s): NP\_058640.1

Human GeneID(s): 50832

Important Comments: The immunizing peptide represents part of an extracellular

domain.

## **Immunogen**

Peptide with sequence C-PELVTTRNNTSFNIS, from the internal region of the protein sequence according to NP\_058640.1.

Please note the peptide is available for sale.

### **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:128000.

**Western blot:** Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

#### **Species Reactivity**

Tested:

Expected from sequence similarity: Human